

N95 Masks Treatable by Battelle CCDS™



3M™ VFLEX 1804



3M™ 1860(S)



3M™ Aura 1870+



3M™ 8205



3M™ 8000



3M™ 8210



3M™ 9010



3M™ Aura 9210+



Gerson 1730



Gerson 2140C



Moldex 1500 Series



Moldex 2200



Kimberly-Clark 46727

Contain Cellulose – Do Not Treat



ProTech NON27501

Check for Cellulose Acetate in undefined masks by placing filter, shell, and coverweb (not straps or nose foam) in Butyric Acid.

Polypropylene and Polyester have satisfactory chemical resistance to it, but Cellulose Acetate (CA) does not so degradation should be evident if the mask contains CA.